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 **Dr. JUNGHANS**[®]
MEDICAL-PRODUKTE • BABYARTIKEL

ENG

Sock and stocking aid of terry cloth for surgical or support stockings

Purpose:

This sock and stocking aid is designed to facilitate the putting on of surgical or support stockings.

Preparation:

- Grip your stockings up to the toes and insert your thumbs in the stocking.
- Now, pull the stocking over the compressed aid.
- Here, make sure the toes are tight over the narrow end of the stocking aid. The remaining stocking should be in the recesses at the upper end.

Use:

- Sit down and allow the prepared stocking aid with the ribbons to slide down onto the floor such that the opening faces you. The stocking aid should rest on the round end. Keep the ribbons in your hands.
- Slide your toes into the stocking aid by moving the ball of your foot forward bit by bit. Be sure to put a little pressure into the downward movement in the process.
- Once your toes are in the tip of the stocking, use the ribbons to pull the stocking aid up via your heel.

Please note:

Apply a little bit of talcum powder for better movement of your foot in the stocking.

Care and cleaning:

Use off-the-shelf cleaning agents to clean the sock and stocking aid.

Guarantee:

Warranty repairs will be available for all defects of this product that are clearly caused by faults in workmanship or material. The warranty period will be twenty-four (24) months. Defects caused by natural wear, intention, negligent or improper handling or use will be excluded.

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