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1 | Introduction

Thank you for deciding for this massager.

This appliance has been tested according to CE and complies with the recognised regulations for this technology and the corresponding safety requirements. This is a quality product that should last a long time when used properly.

2 | Important safety hints

Retain the manual for future use and forward it to other users.

- Do not use the massager:
 - if you have thrombosis
 - varicose veins
 - inflammation of the veins
 - if you have pains of unknown origin
 - if you feel uncomfortable
 - with acute complaints
 - infected skin areas
 - with high fever
 - with an abnormal blood pressure

Ask your doctor.

- Read the instruction manual carefully.
- This device is not intended to be used unsupervised by children, individuals with sensory deprivation of any kind or mentally challenged individuals. Any users falling into these categories should be assisted by a responsible adult.
- Children should not be let without supervision with the appliance in order to avoid them playing with it.
- Please don't use it in the bathroom or other place with high humidity.

3 | Intended use

It is however not intended to treat medical conditions.

The article was developed for massage of calf and arm muscles. Vibration massage to loosen tension or cramping while simultaneously increasing blood circulation in the foot region.

3.1 | Main features of this device

- Extra-soft plush synthetic
- 3 stage selector
- Appealing design

4 | Components

- 1 Calf Massager
- 2 battery compartment
- 3 Hand selector
- 4 pocket for battery compartment

5 | Operating instructions

- Do not submerge the appliance in water or liquid.
- Children should not be let without supervision with the appliance in order to avoid them playing with it.
- If this device is stored near freezing, allow it to acclimate to room temperature before use.

a Installing or replacing batteries:

Open the zipper in the middle of the plush cover. Pull out the battery compartment. Take care not to pull apart the connecting wire.

b Place four AA batteries in compartment in correct position,.

Take note of the correct polarities.

c Lay the battery compartment inside the product.

d Close zipper to the inside padding and outside cover. The item is now ready for use.

CAUTION!

- To prevent damage from battery leakage, remove them if not used for a longer period of time.
- Do not dispose of batteries in fire.
- Make sure the room temperature is right.

6 | Handling

With the selector one may choose among 3 different vibration speeds.

- Press once to turn the machine on and select the lower intensity level.
- Press once more and the machine reaches the higher intensity level.

- With another push of the selector switch the device moves to the highest speed.
- To turn off the machine, press the switch a fourth time.

7 | Uses

- For a soothing massage of the calves, it is best to place both calves are placed on the device.
- Arms and neck may be massaged as shown in the drawing.
- In the regions of the neck, however, a massage should not exceed approximately 5 minutes.

8 | Cleaning and storing

- Use a soft, slightly damp cloth to wipe.
- The product is not machine-washable.
- When cleaning, do not submerge in liquids.
- Do not use gasoline, window cleaner, glass cleaner, furniture polish or paint thinner.
- Store in a clean, dry place.
- Avoid high temperature and direct sunlight.

9 | Delivery Contents

- Calf Massager
- Instruction manual
- Four AA batteries

9.1 | Technical data

Model:	LLU-05
Input voltage:	6V
Battery type:	4 x AA
Material:	95% Cotton, 5% Elastane
Weight:	350 g
Product dimensions:	45 x 19 x 10 cm

10 | Disposal

Electrical and electronic devices may not be disposed of with domestic waste. Consumers are obliged by law to return electrical and electronic devices at the end of their service lives to the public collecting points set up for this purpose or point of sale. Details to this are defined by the national law of the respective country. This symbol on the product, the instruction manual or the package indicates that a product is subject to these regulations. By recycling, reusing the materials or other forms of utilising old devices, you are making an important contribution to protecting our environment.

11 | Warranty

The guarantee lasts for 24 months from the date of purchase for manufacturers' error in material and quality. If it is used for another purpose or operated incorrectly, no liability for possible damages can be accepted. This guarantee is only valid if the instructions in the "instructions for use" have been followed. Warranty is void if product is misused, or repairs attempted by unqualified persons.

Dear Customer,

in case you product does not function as you expected please do not return immediately. Perhaps our hotline [redacted] can be of assistance. If you are still not satisfied, you have the possibility to return products under 40 € value free of charge within the BRD. We require your address. This you may provide to the above listed telephonic hotline or by email to: [redacted]. Prepare the package for shipping and include the receipt and description of the problem. Please be aware of the fact that collect packages may not be accepted.

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