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**Warning:** Please read this instruction manual carefully before using the instrument and keep it for future use.

#### **A. SAFETY INSTRUCTIONS**

1. The device is intended exclusively for private and not for commercial use. Only use this device as describe in this instruction manual.
2. The unit is for personal use only and is not intended for medical or commercial applications.
3. Please note that technically related measuring tolerances are possible, as the scale is not calibrated for professional medical applications.
4. On delivery, the scale is set to “cm” and “kg” . At the back of the scales is a switch to change between “inches” , “pounds” .
5. Place the scale on a firm level floor; a firm floor covering is required for correct measurement.
6. Repairs may only be carried out by professional.
7. Before submitting any complaint, first check the batteries and replace them if necessary.

#### **B. USER-FRIENDLY TIPS**

This digital diagnostic scale is used to weigh you and to provide diagnostic data about your personal fitness.

The scale is equipped with an innovative cordless display that displays the measurements within a range of up to 2 m. The weight is displayed in 100-gram increments. The scale is intended for individual household use. The scale has the following diagnostic functions that can be used by up to 12 persons:

- Body weight measurement
- Measurement of body fat percentage

P 9350 / EF935

- Body water percentage

This scale also has the following functions:

- Switch between kilograms and pounds (LB).
- Automatic shutoff function
- Battery change indicator for weak batteries
- When switched off (standby mode), the scale displays the time.

The scales must not be used by persons with medical implants (e.g. heart pacemakers) otherwise their function could be impaired.




Do not use during pregnancy.

Do not stand on the outer edge to one side of the scale otherwise it may tip!

Keep children away from packaging materials (risk of suffocation).

**This scale is unsuitable for measuring body water in order to draw medical conclusions concerning age related water retention, for example. If necessary ask your physician. Basically, a high body water percentage should be the aim.**


### C. HOW TO USE IT

Note: 1. Please take out the read-out display when setting time.  
 2. Please do not cover the infrared projector and receiver when weighing under wireless transmission mode.  
 3. On wired read-out display, read-out display should be aligned with infrared projector.

### 1) Insert or replace batteries

- Remove the battery cover on the back of the unit.
- Insert 1 x CR2032 batteries observing the correct polarity indicated.
- Replace the cover.
- Please repeat the operation into the scale and control unit.


### 2) Set the time

From standby mode (time displayed), press the [SET] button for more than 3 seconds. The display flashes.


Set the time by pressing the ▲ or ▼ buttons.

Confirm each entry by pressing then [SET] button.

### 3) Set the weight unit

Step onto the scale. The display "0.0 KG" starts flashing. Press [SET] for more than 3 seconds and step off the scale.

- When setting kg as weight unit, the default height unit is cm.


**D. WEIGHTING**

If possible, always weigh yourself at the same time of day (preferably in the morning), after using the toilet, on an empty stomach and unclothed, in order to achieve comparable results.


Important for the measurement: only measure body fat when barefoot and with the soles of the feet slightly moist. Completely dry soles can result in unsatisfactory measurements due to inadequate conductivity.

**INSTALLING THE CONTROL UNIT**

You can either keep the control unit on the scale or hold the control unit in your hand, or even set it on a level surface.


**2) Hold unit in your hand**

Please note with the following:

1. Please do not cover the infrared receiver by hands; Please do not cover the infrared projector by feet.
2. Pay attention to the angle between read-out display and scale body. (Below illustration for your reference) Please don't keep read-out display horizontally, but keep the infrared receiver downward vertically.

### 3) Wall mounting


Using the mounting kit, you can mount the control unit on a wall. Use the mounting kit provided and install the wall holder as shown.

Drill 2 holes in the wall next to each other horizontally.


Insert the wall plugs.

Turn the fastening screws into the wall plugs, until only about 4 mm of the screw heads protrude.

Fasten the wall holder, then set the control unit on the support holder.

Keep scale close to the wall and get better reception.

### E. WEIGHTING WITH DIAGNOSTIC

#### **ADVICES:**

Important for the measurement: only measure body fat when barefoot and with the soles of the feet slightly moist. Completely dry soles can result in unsatisfactory measurements due to inadequate conductivity.

Stand still and upright during the measurement.

Wait a few hours after unaccustomed physical activity.

Wait approx. 15 minutes after getting up in the morning to allow the water in the body to distribute.

Remember that only long-term trends are important. Short-term changes in weight within a few days are usually caused by a loss of fluid. Body water plays an important role towards our general wellbeing.

#### **LIMITATIONS:**

When measuring body fat and other values, deviating and implausible results may occur in:

- children under approx. 10 years,
- competitive athletes and body builders,
- pregnant women,
- persons with fever, undergoing dialysis treatment or with symptoms of edema or osteoporosis,
- persons taking cardiovascular medication (affecting the heart and vascular system),
- persons taking vascular dilating or vascular constricting medication,
- persons with considerable anatomic deviations of the legs in relation to total body size (length of the legs considerably shortened or lengthened).

#### **BIOELECTRIC IMPEDANCE ANALYSIS PRINCIPLE:**

This scale operates according to the B.I.A. principle (bioelectric impedance analysis). This enables the measurement of physical relationships within seconds by means of an undetectable, completely harmless electric current. The body fat percentage and other physical relationships in the body can be determined by measuring the electrical resistance (impedance) and calculating constants and individual parameters such as age, height, AND gender.

Muscle tissue and water have good electrical conductivity, and therefore lower resistance. Bones and fatty tissue, on the other hand, have low conductivity as fat cells and bones conduct hardly any current as a result of their very high resistance.

Please be aware that values obtained from the diagnostic scale represent only an approximation of actual analytical medical data. Only a specialist physician can accurately determine body fat, body water percentage using medical procedures (e.g. computed tomography).


### 1) Setting user data

You must enter personal user parameters in order to measure your body fat percentage and other physical data.

The scale has 12 user memory locations which allows, for example, you and your family members to save personal settings and retrieve them later.


Important:

Note that there may be no contact between your two feet, legs, calves or thighs. Otherwise the measurement will not be performed correctly.

User data	Parameters
Memory position	1~12
Body height	100 ~250 cm (3'03.5" ~8'02.5")
Birthday	year, month, day
Sex	male ♀, female ♂

### 2) Weighting

The scale always switches off automatically after a certain time. The time appears in the display.

**3) Reading results:**

Age	Fat(F)	Hydration(F)	Fat(M)	Hydration(M)	Status	Key Line Indication
≤30	4.0-16.0%	66.0-57.8%	4.0-11.0%	66.0-61.2%	Too lean	
	16.1-20.5%	57.7-54.7%	11.1-15.5%	61.1-58.1%	Lean	
	20.6-25.0%	54.6-51.6%	15.6-20.0%	58.0-55.0%	Normal	
	25.1-30.5%	51.5-47.8%	20.1-24.5%	54.9-51.9%	Fat	
	30.6-60.0%	47.7-27.5%	24.6-60.0%	51.8-27.5%	Too fat	
>30	4.0-20.0%	66.0-55.0%	4.0-15.0%	66.0-58.4%	Too lean	
	20.1-25.0%	54.9-51.6%	15.1-19.5%	58.3-55.3%	Lean	
	25.1-30.0%	51.5-48.1%	19.6-24.0%	55.2-52.3%	Normal	
	30.1-35.0%	48.0-44.7%	24.1-28.5%	52.2-49.2%	Fat	
	35.1-60.0%	44.6-27.5%	28.6-60.0%	49.1-27.5%	Too fat	

• The above information is only for reference


**F. TROUBLE SHOOTING**

Low battery power for read-out display. Please replace battery.

Low battery power for main unit. Please replace battery.

Over-load Indication  
The weighing subject on the platform exceeds the maximum capability of scale. Please step off to avoid damage.

Low Fat% Indication  
The Fat% is too low. Please lead a more nutritional diet and take good care of yourself.

High Fat% Indication  
The Fat% is too high. Please watch your diet and do more exercise.

Transmission error. Check that your feet or hand do not interfere with the infrared transmission and that the receiver is facing down and towards the transmitter.

**G. CARE AND MAINTENANCE**

The precision of the measured values and service life of the unit depend on careful use. The unit should be cleaned occasionally. Do not use abrasive detergents and never immerse the unit in water.

Ensure that no liquid gets on the scale. Never dip the scale into water. Never wash it under running water.

Clean the product only with a dry or slightly damp, lint-free cloth.

Do not use abrasive cleaners to clean the unit.

Do not place any objects on the scale when it is not being used.

Protect the scale from impacts, moisture, dust, chemicals, large changes in temperature and sources of heat that are too close (oven, radiator).

Do not press the buttons using excessive force or with pointed objects.

Do not expose the scale to high temperatures or strong electromagnetic fields (e.g. mobile phones).

Remove batteries if the unit will not be used for an extended period of time. This prevents damage that may occur from leaks.

**H. FEATURES**

Material: Tempered glass top, ABS and HIPS

Power supply: 2 x 3V CR 2032 batteries (included)

Infrared Technology

Maximum load: 150 kg - Graduation: 100 grams

Screen size: 46.5 x 27.5 mm

The European directive 2012/19/EC about electronic and electric waste, requires that you can't throw away defective domestic appliance with common waste. Used device has to be collect separately to optimize the recuperation rate and the recycling of the materials to reduce the impact on the health and the environment.

This appliance is not intended for use by persons (including children) with reduced physical, sensory or mental capacities, or by persons lacking in experience or knowledge, unless they are supervised or instructed on the use of this apparatus by a person responsible for their safety. Children should be supervised to ensure they do not play with this appliance.

**Precision and precautions about the battery use**

- Throwing away a used battery in the environment or with household waste pollutes and avoids the recuperation of reusable materials.
- It is important to limit the battery consumption and to respect the following rules:
  - Favor alkaline battery (last longer than saline battery) and when it's possible use rechargeable battery.
  - Leave used batteries and accumulators in adapted containers that you can find in some retailers. In that way, metals don't pollute the environment and the health.

P 9350 / EF935

- *Battery has to be inserted by respecting the polarity show on the apparel and the battery. An incorrect positioning can damage the apparel, cause leaking and even provoke fire and the explosion of the battery.*
- *To assure a good functioning, battery has to be in a good condition. If there is a problem, put new battery.*
- *Never try to recharge battery not made for this. It could leak, warm up, cause fire or explode.*
- *Replace all the battery at the same time. Never mix up alkaline battery with saline battery.*
- *Used battery has to be removed.*
- *Also, remove the battery if you're not using the apparel for a long time.*
- *Never try to short-circuit the terminal of a battery.*
- *Never throw away the battery in fire because they can explode.*
- *Take out the accumulator from the apparel before reload it.*
- *We recommend that an adult watch if a child change the battery to respect those rules.*
- *If you swallow a battery, go immediately to your doctor or the hospital.*
- *For the best ongoing performance, we recommend using only fresh, high-quality, alkaline batteries.*

**Importé par / Importiert von / Geïmporteerd door / Imported by :**  
**CAPTELEC, 59170 Croix - FRANCE**

Carton et emballage papier à  
trier / Kartons und  
papierverpackung recyceln /  
Verpakkingsmateriaal in karton  
en papier te sorteren / Recycle  
carton and paper packaging